

INDICATOR OF REMAINING ENERGY IN STORAGE CELL OF IMPLANTABLE MEDICAL DEVICE

ABSTRACT

5 A manganese dioxide (MnO₂) or silver vanadium oxide (SVO) or other
battery of an implantable medical device having a relatively flat quiescent battery
voltage during a beginning portion of the battery's useful life, makes it difficult to
use quiescent battery voltage as an indicator of remaining battery energy during this
portion of the battery life. A substantially constant load current pulse is drawn from
10 the battery and a pair of loaded battery terminal voltage measurements is taken
during this pulse. A difference between the voltage measurements is computed.
This difference can be expressed as a rate of change, a slope, or a polarization angle,
and can be used with stored data from similar batteries to determine remaining
energy of the battery. A quiescent battery voltage can also be used in combination
15 with this technique, and/or for distinguishing between different remaining energies
corresponding to the same difference, slope, or polarization angle.